



The only automated
chest compression system
improved specifically with
the cath lab in mind

New **radiotranslucent**
back plate



Now, with a new carbon fiber back plate, visualizing coronary arteries during ongoing compressions is **even easier.**



The LUCAS cath lab back plate is made out of carbon fiber which is lightweight, very strong, and fully radiotranslucent. The material enables the construction of a thin but durable back plate that does not interfere with the angiogram. LUCAS facilitates the management of acute and life-threatening complications in the cath lab:

- **Designed to deliver effective, uninterrupted and consistent compressions for victims of circulatory arrest, according to AHA guidelines**
- **Facilitates circulation:**
 - to vital organs
 - by allowing compressions on a patient on the cath lab table
 - by minimizing interruptions despite ongoing fluoroscopy
- **Allows for simultaneous catheterization as well as PCI during ongoing CPR**
- **Frees up resources and buys time for decision making in a stressful situation**
 - a bridge to continued PCI, LVAD, cardiopulmonary bypass or surgery
 - increases personnel safety by minimizing radiation exposure typically seen during manual CPR



“The mood in the cath lab was calm at all times despite the ongoing VF... This is quite contrary to what usually happens in such situations when manual compressions are used.”

– Olivecrona, Lund, Sweden (tctmd.com 24 Oct 2006).

“We believe LUCAS CPR is a big advancement in resuscitation therapy. Furthermore, it enables for a fast intervention of a closed coronary artery also during instable circulation, and allows for a shortened door-to-balloon time.”

– Schäfer et al, Dresden Germany (*Clin Res Cardiol* 96: Suppl 1 2007)

LUCAS and Fluoroscopy



LUCAS is radiotranslucent, except for the hood and the piston. This allows for these projections in monoplane

**Straight
Cranial**



**Straight
Caudal**



**Straight
Lateral**



**LAO/RAO
Cranial**



**LAO/RAO
Caudal**



PCI Back Plate Specifications

- Material: Carbon Fiber
- Dimensions: (L x W) : 49 x 24 cm / 19.3 x 9.4 inches
Thickness : 7 mm / 0.3 inches
- Weight: 1 kg / 2.2 lbs
- Fully radiotranslucent
- Can be pre-positioned in unstable patients, without compromising angiographic imaging, and allowing for quick set up in case of refractory/non-shockable arrest.

”... percutaneous intervention can be carried out successfully during ongoing automatic mechanical compressions by the LUCAS device, which provides an efficient circulatory pulsatile support [...] without the need for additional staff involved in basic life support.

– Agostoni et al, Antwerp, Belgium (*Int J Cardiol.* 2007 Feb 28).

LUCAS is manufactured by JOLIFE AB in Sweden and distributed worldwide by Physio-Control, Inc.

For further information, please contact Physio-Control at 800.442.1142 (U.S.), 888.879.0977 (Canada) or visit our website at www.physio-control.com.



Physio-Control Headquarters

11811 Willows Road NE
Redmond, WA 98052
www.physio-control.com

Customer Support

P. O. Box 97006
Redmond, WA 98073
Toll Free 800 442 1142
Fax 800 426 8049

Physio-Control Canada

Medtronic of Canada Ltd
99 Hereford Street
Brampton, ON
L6Y 0R3
Tel 888 879 0977
Fax 866 430 6115



Jolife AB, Scheelevägen 17, SE-223 70 LUND, Sweden